

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

BLUE BHD AUTOCLAVE BAG WITH TYVEK® DRAWSTRING



TYPICAL APPLICATION & CHARACTERISTICS:

- Designed for steam sterilization
- Incremental improvement compared to cutting and wrapping with Bio-Shield™ and autoclave tape
- Avoid sterilization Re-validation
- Tyvek® drawstring closure
- Steam sterilization indicator
- Eliminates need to cut material and the use of particulate generating autoclave tape
- Fast and consistent
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Resists tearing and punctures
- 2-ply for good microbial barrier
- Manufactured and packaged in ISO certified cleanrooms
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	BHD MATERIAL
Basis Weight	Gm/m2	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters/ 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams/sq. cm	2,150
Flammability	seconds	8

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

BLUE BHD AUTOCLAVE BAG WITH TYVEK® DRAWSTRING



KEYSTONE PART NUMBER	DIMENSIONS
BNWBG0642D-2P	6" X 42"
BNWBG1212D-2P	12" X 12"
BNWBG1624D-2P	16" X 24"
BNWBG1644D-2P	16" X 44"
BNWBG1830D-2P	18" X 30"
BNWBG1854D-2P	18" X 54"
BNWBG2436D-2P	24" X 36"
BNWBG4648D-2P	46" X 48"