

SINCE 1925

# KEYSTONE

## SAFETY

# FACE MASK WITH EAR LOOP

## DESCRIPTION

Face masks with ear loops and adjustable nose clip securely covers the nose and mouth. Disposable, single-use face masks are hypoallergenic, latex-free and non-sterile. Its pleated, non-woven material structure has a full width adjustable nose piece provides a Bacterial Filtration Efficiency (BFE) of 95% or greater.

## CHARACTERISTICS

- 100% Latex-free
- 3 Layer construction
- Hypoallergenic
- Pleated Design

## FM-EL-BLUE



**SIZES AVAILABLE:** One size fits all

**COLORS AVAILABLE:** White, Blue, Pink

**PACKAGING:** 50/dispenser box-40 bxs/cs = 2000

**COUNTRY OF ORIGIN:** China

All Keystone Safety Product components comply with the provisions of the Federal Food, Drug and Cosmetic Act, which allow repeated use in direct contact with food. Fabric conforms to FDA code of Federal Regulations Title 21, Rule 110.10 & 177.1950

ITEM	SIZE	COUNT	WT/BOX	DIMENSION	CUBE	QUANTITY/ SKID	SKID CONFIG TIERS x ROWS
FM-EL-BLUE	OSFA	2000	17.1	21x16x12	.062	48	6x8

1591 HYLTON ROAD PENNSAUKEN, NJ 08110  
KEYSTONESAFETY.COM | INFO@KEYSTONESAFETY.COM

TOLL FREE (800) 663.5439  
TOLL FREE FAX (800) 663.5432  
PHONE (856) 663.5740  
FAX (856) 663.6075

SINCE 1925

# KEYSTONE

## SAFETY