


SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

SMS STERILIZATION WRAPPER HEAVY DUTY 50 GSM



Keystone Cleanroom Products' line of Heavy Duty SMS wrappers provide an excellent microbial barrier. Ideal for pharmaceutical production applications.

- Manufactured and packaged in a certified ISO8 Cleanroom
- Clean--no cellulose
- Strong, rip and tear resistant 50 gsm SMS material
- Soft and drape-able fabric
- An extensive range of sizes available means less cutting
- Hydrophobic barrier allows sterilization gasses, but prevents liquid and microbial penetration
- Sub packed into sealed labeled poly bags with Lot Number and Date of Manufacture for complete traceability
- Certificate of Conformance sent with each shipment
- Made in the USA

For use in steam, ethylene oxide and hydrogen peroxide sterilization processes.

MATERIAL SPECIFICATIONS:


PROPERTY	SPECIFICATION	UNIT
Basis Weight	50	Gm/m ²
Grab Strength (CD)	20	lbs
Tear Strength (MD)	7.5	lbs
Mullen Burst Strength (Gross)	34	psi
Air Permeability	50	cfm
Hydrohead	55	mbar

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

SMS STERILIZATION WRAPPER HEAVY DUTY 50 GSM



KEYSTONE PART NUMBER	DIMENSIONS	KEYSTONE PART NUMBER	DIMENSIONS
KCRP-300-01.5X01.5	1.5" X 1.5"	KCRP-300-2020	20" X 20"
KCRP-300-0202	2" X 2"	KCRP-300-2424	24" X 24"
KCRP-300-0303	3" X 3"	KCRP-300-3030	30" X 30"
KCRP-300-0404	4" X 4"	KCRP-300-3636	36" X 36"
KCRP-300-0606	6" X 6"	KCRP-300-4040	40" X 40"
KCRP-300-0808	8" X 8"	KCRP-300-4545	45" X 45"
KCRP-300-0909	9" X 9"	KCRP-300-4548	45" X 48"
KCRP-300-0912	9" X 12"	KCRP-300-4848	48" X 48"
KCRP-300-1212	12" X 12"	KCRP-300-5454	54" X 54"
KCRP-300-1515	15" X 15"	KCRP-300-5472	54" X 72"
KCRP-300-1818	18" X 18"	KCRP-300-6060	60" X 60"
KCRP-300-1822	18" X 22"	KCRP-300-6080	60" X 80"