

SINCE 1925

KEYSTONE

SAFETY

CPE ISOLATION GOWN

DESCRIPTION

CPE (Cross Linked Polyethylene) isolation gowns are FDA compliant and offer excellent barrier protection in many applications. These gowns have long sleeves, a rounded pull over neck, ties at waist and thumb loop cuffs. CPE gowns are made from non-toxic, food safe, waterproof material, which is commonly used in the food manufacturing and processing, industrial factories, medical and hospital facilities. It is ideal for economical protection against water, dust and splashes.

CHARACTERISTICS

- 100% latex free
- Produced in an ISO 9001 certified factory

AAMI LEVEL 2 PROTECTION

This isolation gown meets the Association for the Advancement of Medical Instrumentation (AAMI) PB70:2012 standard. "Liquid barrier performance and classification of protective apparel and drapes intended for use in healthcare facilities."

SIZES AVAILABLE: 47" and 55"

COLORS AVAILABLE: White, Blue, Yellow, Green

PACKAGING: 100/case

COUNTRY OF ORIGIN: China

All Keystone Safety Product components comply with the provisions of the Federal Food, Drug and Cosmetic Act, which allow repeated use in direct contact with food. Fabric conforms to FDA code of Federal Regulations Title 21, Rule 110.10 & 177.1950

ISO-TL



| ITEM | SIZE | COUNT | WT/BOX | DIMENSION | CUBE | QUANTITY/ SKID | SKID CONFIG TIERS x ROWS |
|--------|------|-------|--------|-----------|------|-------------------|-----------------------------|
| ISO-TL | 47" | 100 | 10.6 | 14x9x15 | .029 | 54 | 9x6 |
| ISO-TL | 55" | 100 | 12.2 | 14x9x15 | .029 | 54 | 9x6 |

1591 HYLTON ROAD PENNSAUKEN, NJ 08110
KEYSTONESAFETY.COM | INFO@KEYSTONESAFETY.COM

TOLL FREE (800) 663.5439
TOLL FREE FAX (800) 663.5432
PHONE (856) 663.5740
FAX (856) 663.6075

SINCE 1925

KEYSTONE

SAFETY