

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PRE-STERILIZED AUTOCLAVABLE BIOHAZARD BAGS



Keystone Cleanroom Products Gamma Irradiated Biohazard Bags are useful to decontaminate waste or potentially hazardous materials that are exiting the aseptic production area through the autoclave.

The irradiated biohazard bags eliminate the introduction of contamination into the aseptic production area..

TYPICAL APPLICATION & CHARACTERISTICS:

- Packaged in an ISO6 certified Cleanroom
- Strong, 2 mil clear polypropylene
- Validated Sterile to a 10^{-6} SAL
- **Autoclavable up to 135° C (275° F)**
- Certificates of Conformance and Sterility sent with each shipment
- Conveniently packaged in Tyvek®/Mylar chevron peel pouches: 1, 2, 5 or 10 per pack configurations
- Made in USA

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PRE-STERILIZED AUTOCLAVABLE BIOHAZARD BAGS

Keystone CRP Part No.	Description
1015CBHNPST1	10" x 15", 1 per pack
1015CBHNPST2	10" x 15", 2 per pack
1015CBHNPST5	10" x 15", 5 per pack
1015CBHNPST10	10" x 15", 10 per pack
1224CBHNPST1	12" x 24", 1 per pack
1224CBHNPST2	12" x 24", 2 per pack
1224CBHNPST5	12" x 24", 5 per pack
1224CBHNPST10	12" x 24", 10 per pack
1924CBHNPST1	19" x 24", 1 per pack
1924CBHNPST2	19" x 24", 2 per pack
1924CBHNPST5	19" x 24", 5 per pack
1924CBHNPST10	19" x 24", 10 per pack
2430CBHNPST1	24" x 30", 1 per pack
2430CBHNPST2	24" x 30", 2 per pack
2430CBHNPST5	24" x 30", 5 per pack
2430CBHNPST10	24" x 30", 10 per pack
2436CBHNPST1	24" x 36", 1 per pack
2436CBHNPST2	24" x 36", 2 per pack
2436CBHNPST5	24" x 36", 5 per pack
2436CBHNPST10	24" x 36", 10 per pack
3138CBHNPST1	31" x 38", 1 per pack
3138CBHNPST2	31" x 38", 2 per pack
3138CBHNPST5	31" x 38", 5 per pack
3138CBHNPST10	31" x 38", 10 per pack
3748CBHNPST1	37" x 48", 1 per pack
3748CBHNPST2	37" x 48", 2 per pack
3748CBHNPST5	37" x 48", 5 per pack
3748CBHNPST10	37" x 48", 10 per pack