

SINCE 1925

**KEYSTONE**

SAFETY ↔ CLEANROOM PRODUCTS

# PRE-STERILIZED MOP BUCKET LINERS



The use of Gamma Irradiated Mop Bucket Liners eliminates the need to clean and sterilize mop buckets before and after each use. The liners will extend the lifespan of stainless-steel buckets. Clean on your own schedule without having to wait for janitorial autoclave loads to be completed. This frees up valuable autoclave resources and reduces the utility costs on expensive steam and WIFI. The use of the liners minimizes the risks of spreading mold spores during cleaning.

## TYPICAL APPLICATION & CHARACTERISTICS:

- Packaged in ISO certified Cleanroom
- Strong 4-mil polyethylene
- Standard sizes are 28" x 28" and 34" x 28"
- Conveniently packaged in Tyvek®/Mylar chevron peel pouches; 1, 2, 3, 4 per pack configuration or
- **Validated Sterile to a 10<sup>-6</sup> SAL**
- **Certificates of Conformance and Sterility sent with all shipments**
- Made in USA

## MATERIAL SPECIFICATIONS:

PHYSICAL PROPERTIES	UNITS
Material-Low Density Polyethylene	4 mil
Dimensions	(+/- 1")
S-2828CP4	28" X 28"
S-3428CP4	34" X 28"

SINCE 1925

**KEYSTONE**

SAFETY ↔ CLEANROOM PRODUCTS

# PRE-STERILIZED MOP BUCKET LINERS



KEYSTONE PART NUMBER	DIMENSIONS
S-2828CP4-1	28" X 28", 1 per pack
S-2828CP4-2	28" X 28" , 2 per pack
S-2828CP4-3	28" X 28" , 3 per pack
S-2828CP4-4	28" X 28" , 4 per pack
S-3428CP4-1	34" X 28" , 1 per pack
S-3428CP4-2	34" X 28" , 2 per pack
S-3428CP4-3	34" X 28" , 3 per pack
S-3428CP4-4	34" X 28" , 4 per pack