


SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

ONE-STEP, 2-PLY SUPER DUTY SMS STERILIZATION WRAP



Keystone Cleanroom Products' line of One-Step, Super Duty SMS 2-ply wrappers provide an excellent microbial barrier. Ideal for pharmaceutical production applications where two layers are required. Saves time and reduces the amount of autoclave tape.

- Manufactured and packaged in a certified ISO6 Cleanroom
- Clean--no cellulose
- Strong, rip and tear resistant 88 gsm SMS material
- Soft and drape-able fabric
- An extensive range of sizes available means less cutting
- Hydrophobic barrier allows sterilization gasses, but prevents liquid and microbial penetration
- Sub packed into sealed labeled poly bags with Lot Number and Date of Manufacture for complete traceability
- Certificate of Conformance sent with each shipment
- Made in the USA

For use in steam, ethylene oxide and hydrogen peroxide sterilization processes.

MATERIAL SPECIFICATIONS:

PROPERTY	SPECIFICATION	UNIT
Basis Weight	88	Gm/m2
Grab Strength (CD)	40	lbs
Tear Strength (MD)	11	lbs
Mullen Burst Strength (Gross)	60	psi
Air Permeability	25	cfm
Hydrohead	70	cm

*The information above references the specifications of each layer of the material. The One-Step Sterilization wrap is made with two plies of material, ultrasonically bonded on opposite edges.

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

ONE-STEP, 2-PLY SUPER DUTY SMS STERILIZATION WRAP



KEYSTONE PART NUMBER	DIMENSIONS	KEYSTONE PART NUMBER	DIMENSIONS
KCRP1-600-0808	8" x 8" 20 cm x 20 cm	KCRP1-600-3636	36" x 36" 91 cm x 91 cm
KCRP1-600-0909	9" x 9" 23 cm x 23 cm	KCRP1-600-4040	40" x 40" 102 cm x 102 cm
KCRP1-600-0912	9" x 12" 23 cm x 30 cm	KCRP1-600-4545	45" x 45" 114 cm x 114 cm
KCRP1-600-1212	12" x 12" 30 cm x 30 cm	KCRP1-600-4548	45" x 48" 114 cm x 122 cm
KCRP1-600-1515	15" x 15" 38 cm x 38 cm	KCRP1-600-4848	48" x 48" 122 cm x 122 cm
KCRP1-600-1818	18" x 18" 46 cm x 46 cm	KCRP1-600-5454	54" x 54" 137 cm x 137 cm
KCRP1-600-1822	18" x 22" 46 cm x 56 cm	KCRP1-600-5472	54" x 72" 137 cm x 183 cm
KCRP1-600-2020	20" x 20" 51 cm x 51 cm	KCRP1-600-5490	54" x 90" 137 cm x 229 cm
KCRP1-600-2424	24" x 24" 61 cm x 61 cm	KCRP1-600-6060	60" x 60" 152 cm x 152 cm
KCRP1-600-3030	30" x 30" 76 cm x 76 cm	KCRP1-600-6080	60" x 80" 152 cm x 203 cm