

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

TYVEK® PRODUCT TRANSFER SLEEVES



Tyvek® Product Transfer Sleeves move powder between pharmaceutical processing equipment and containers during production. These sleeves can be designed with an “inner sleeve” and a breathable “outer skirt” to ensure that the product is moved precisely into the receptacle, while also minimizing the “dusting” or “blowback” that can occur. These sleeves are disposable, but are rugged enough to be used as dedicated sleeves for a particular product and re-used many times before discarding.

TYPICAL APPLICATION & CHARACTERISTICS:

- Individually poly bagged and labeled for easy identification
- Traceable via manufacturer lot number
- Simple to install by one person, no tools needed
- Breathable Polypropylene top allows for outgassing
- Minimizes waste by product being unnecessarily dispersed throughout the facility
- Reduces operator exposure and room cleaning
- Optional elastic, ties or adjustable rubber roping with a cinch closure, for the most secure fit
- Anti-static Tyvek® material reduces sparking and explosion risks
- Eliminates the need for laundering, cleaning validation and repairs common with cloth sleeves
- Tough Tyvek® fabric resists tearing and punctures
- Lint-free
- Manufactured and packaged in ISO certified Cleanroom
- Certificates of Conformance sent with shipment
- Made in USA

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

TYVEK® PRODUCT TRANSFER SLEEVES



MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A TYVEK®
Basis Weight	g/m2	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M2/24 HR	ASTM 398	1700
Hydrostatic Head	cm H2O	ISO 811	120
Trapezoidal Tear MC/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50cm2)	108
Puncture Resistance	N	EN 863	10.8
Melting Point	C °		135

KEYSTONE PART NUMBER	DESCRIPTION
FITZTWPPBL30	1422A Tyvek® Fitzmill Transfer Sleeve to 30 Gallon Drum, with Polypropylene Breather Panel Top and with Barrel Lock closures, 1 per poly bag
FITZWPP30	1422A Tyvek® Fitzmill Transfer Sleeve to 30 Gallon Drum, with Polypropylene Breather Panel Top, 1 per poly bag
FITZTWPPBL	1422A Tyvek® Universal Fitzmill Transfer Sleeve to 55 gallon Drum, with Polypropylene Breather Panel Top and with Barrel Lock closures, 1 per poly bag
FITZTWBL	1422A Tyvek® Universal Fitzmill Transfer Sleeve to 55 Gallon Drum, with Barrel Lock closures, 1 per poy bag
FITZTWPPNC	1422A Tyvek® Universal Fitzmill Transfer Sleeve to 55 Gallon Drum, with Polypropylene Breather Panel Top, No Clamps, 1 per poly bag
FITZTWNC	1422A Tyvek® Universal Fitzmill Transfer Sleeve to 55 Gallon Drum, No Clamps, 1 per poly bag
FITZTWWWBL	1422A Tyvek® Fitzmill Transfer Sleeve to 55 Gallon Drum with 7" x 12" Window, Barrel Lock closures, 1 per poly bag
FITZTWWWNC	1422A Tyvek® Fitzmill Transfer Sleeve to 55 Gallon Drum with 7" x 12" Window, No Clamps, 1 per poly bag